

DOCUMENT RESUME

ED 153 868

SE 024 252

AUTHOR Bonte, John I.; Davidson, Arnold C.
TITLE Alkalinity Analysis. Training Module 5.220.2.77.
SPONS AGENCY Department of Labor, Washington, D.C.; Iowa State
Dept. of Environmental Quality, Des Moines.
PUB DATE Sep 77
NOTE 33p.; For related documents, see SE 024 249-254;
Contains marginal legibility in transparency
masters

EDRS PRICE MF-\$0.83 HC-\$2.06 Plus Postage.
DESCRIPTORS *Chemical Analysis; Chemistry; *Instructional
Materials; *Laboratory Procedures; *Laboratory
Techniques; *Post Secondary Education; Secondary
Education; Units of Study
IDENTIFIERS *Alkalinity; *Water Treatment

ABSTRACT

This document is an instructional module package prepared in objective form for use by an instructor familiar with the acid-base titrimetric procedure for determining the hydroxide, carbonate and bicarbonate alkalinity of a water sample. Included are objectives, an instructor guide, student handouts and transparency masters. A video tape is also available from the author. This module considers use of the pH meter, preparation and standardization of reagents, titration and calculation of results. (Author/RE)

* Reproductions supplied by EDRS are the best that can be made *
* from the original document. *

ED153868

U.S. DEPARTMENT OF HEALTH,
EDUCATION & WELFARE
NATIONAL INSTITUTE OF
EDUCATION

THIS DOCUMENT HAS BEEN REPRODUCED EXACTLY AS RECEIVED FROM THE PERSON OR ORGANIZATION ORIGINATING IT. POINTS OF VIEW OR OPINIONS STATED DO NOT NECESSARILY REPRESENT OFFICIAL NATIONAL INSTITUTE OF EDUCATION POSITION OR POLICY.

ALKALINITY ANALYSIS
Training Module 5.220.2.77

"PERMISSION TO REPRODUCE THIS MATERIAL HAS BEEN GRANTED BY

Mary Jo Bruett

TO THE EDUCATIONAL RESOURCES INFORMATION CENTER (ERIC) AND USERS OF THE ERIC SYSTEM."

Prepared for the
Iowa Department of Environmental Quality
Wallace State Office Building
Des Moines, Iowa 50319

by

John L. Bonte
Developer
Arnold C. Davidson
Project Director
Clinton Community College
1000 Lincoln Boulevard
Clinton, Iowa 52732

The publication of these training materials was financially aided through a contract between the Iowa Department of Environmental Quality and the Office of Planning and Programming, using funds available under the Comprehensive Employment and Training Act of 1973. However, the opinions expressed herein do not necessarily reflect the position or policy of the U. S. Department of Labor, and no official endorsement by the U. S. Department of Labor should be inferred.

September, 1977

1/2 2

FE 024 252

Module No:	Module Title: Alkalinity Analysis
Approx. Time: 3 hours	Submodule Title: Topic: Summary

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Operate and standardize a pH meter.
2. Determine hydroxide, carbonate, and bicarbonate alkalinity of a water sample.
3. Prepare and standardize all reagents needed for an alkalinity analysis.

Instructional Aids:

Transparencies Ak1 - Ak6 - softening videotape

Instructional Approach:

Lecture, discussion, videotape viewing, laboratory practice.

References:

1. "Standard Methods for the examination of water and wastewater," 14th ed., pp. 278-282.
2. G.V. James, "Water Treatment," 3rd ed., London: Technical Press, 1965.

Class Assignments:

Module No:	Module Title: Alkalinity Analysis
Approx. Time: 0.5 hours	Submodule Title:
	Topic: Definition of Alkalinity

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Write chemical reactions which explain alkalinity.
2. Explain how alkalinity is used in water treatment.
3. Describe the difference between phenolphthalein and total alkalinity.
4. Describe the difference between hydroxide, bicarbonate, and carbonate alkalinity.
5. Relate pH to hydroxide and hydrogen ion concentrations.

Instructional Aids:

Transparency Ak1 - titration end points.

Transparency Ak2 - Relation of pH to concentration.

Softening video tape

Instructional Approach:

Lecture - discussion

References:

1. Standard Methods, p. 278
2. James, p. 107-122.

Class Assignments:

Module No: °	Module Title: Alkalinity Analysis
	Submodule Title:
Approx. Time: 0.25 hours	Topic: Use of the pH meters

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Describe the function of the pH meter.
2. Properly standardize the pH meter.
3. Determine the pH of a solution with a pH meter.

Instructional Aids:

Transparency Ak3 - diagram of pH meter.

Instructional Approach:

Lecture - demonstration

References:

Operator's manual for pH meter used.

Class Assignments:

Module No:
Ak

Topic:
Use of the pH meter

Instructor Notes:

Instructor Outline:

Transparency Ak-3
Diagram of pH meter

1. pH meter measure the H^+ ion concentration of a solution electronically.
2. Standardize with pH 7 buffer solution to a meter reading of 7.0.
3. Determine pH by placing electrodes in solution and reading meter.

Module No:	Module Title: Alkalinity Analysis
Approx. Time: 0.25 Hours	Submodule Title: Topic: Safety

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Locate the following in the laboratory and demonstrate proper use: emergency shower, fire extinguisher, eye wash, first aid kit.
2. Select and use safety glasses, lab coat or apron and gloves in the appropriate situation.
3. Describe the hazards associated with the chemicals used in the alkalinity determination.

Instructional Aids:

Hand-out of safety rules for the laboratory.

Instructional Approach:

Lecture/demonstration

References:

Basic laboratory skills module

Class Assignments:

Module No:
Ak

Topic:
Safety

Instructor Notes:

Instructor Outline:

1. Point out to student all the safety equipment in the laboratory.
2. Point out the hazard of strong acids in the eye.
3. Electrical shock from pH meter and acid burns are the main hazards.

Module No:	Module Title: Alkalinity Analysis
	Submodule Title:
Approx. Time: 0.5 hours	Topic: Preparation of reagents

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Prepare a 0.05 standard sodium carbonate solution and calculate its exact normality.
2. Prepare 0.1N and 0.02N hydrochloric acid solutions.

Instructional Aids:**Instructional Approach:**

Laboratory practice

References:

Standard Methods, p. 279

Class Assignments:

Module No: Ak	Topic: Preparation of reagents
Instructor Notes:	Instructor Outline:
<p>1. 2.5g per liter</p> <p>2.a. 0.1 N HCl 8.3 ml conHCl per liter</p> <p>b. 0.02 N HCl 200 ml 0.1 NHCl per liter of solution</p>	<p>1. Na_2CO_3 solution: N. =0.05 N $N = \frac{g \text{ Na}_2\text{CO}_3}{53} \text{ per liter}$</p> <p>2. Acid a. 0.1N hydrochloric acid b. 0.02 N hydrochloric acid</p>

Module No:	Module Title: Alkalinity Analysis
Approx. Time: 0.5 hours	Submodule Title: Topic: Standardization of Acid

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Properly titrate the standard Na_2CO_3 solution with 0.1N HCl potentiometrically and construct a titration curve.
2. Identify inflection points.
3. Calculate from the titration curve the exact normality of the 0.1N HCl and the 0.02N HCl.
4. Calculate the CaCO_3 equivalence of the HCl solutions.

Instructional Aids:

Transparency Ak4 - sample Na_2CO_3 titration curve.

Instructional Approach:

Laboratory practice

References:

Standard Methods p. 279

Class Assignments:

Module No:

Ak

Topic:

Standardization of acid

Instructor Notes:

Instructor Outline:

Transparency Ak-4
 Na_2CO_3 titration curve

1. Have students measure ml acid (.1N) vs. pH and plot points. Titrate 40 ml Na_2CO_3 solution.
2. Point out inflection points on standard curve.
3.
$$N(.1\text{NHCl}) = \frac{\text{ml Na}_2\text{CO}_3}{\text{ml HCl}} \times (N \text{ Na}_2\text{CO}_3)$$

$$N(.02 \text{ N HCl}) = N(.1 \text{ NHCl}) \times 0.2$$
4. CaCO_3 equivalence (mg/ml) = 50 X N(HCl)

Module No:	Module Title: Alkalinity Analysis
Approx. Time: 0.5 hours	Submodule Title:
	Topic: Titration of water samples

Instructional Objective:

Upon completion of this module the participant should be able to:

1. Properly titrate a water sample potentiometrically with standard HCl, construct a titration curve, and identify inflection points on the curve.
2. Properly titrate a water sample to a predetermined pH and record appropriate data.
3. Titrate a water sample of low alkalinity to determine total alkalinity.

Instructional Aids:

Transparency Ak5 - sample titration curve

Instructional Approach:

Laboratory practice

References:

Standard Methods, p. 280.

Class Assignments:

Module No: Ak	Topic: Titration of water samples
Instructor Notes:	Instructor Outline:
Transparency Ak-5 Sample titration curve	<ol style="list-style-type: none">1. Students perform titration and plot points. Point out inflections on sample curve.2. Alternatively have students titrate to predetermined pH(8.3 and 3.7-5.1) and compare with inflection method.3. For low alkalinity:<ol style="list-style-type: none">a. Use .02 N acid, 100 ml sampleb. titrate to 4.5 - record exact pH and volume added(B ml)c. titrate to exactly 0.3 pH units lower, record volume(C ml) <p>Total Alkalinity as CaCO_3=</p> $\frac{(2B-C) \times N(\text{HCl}) \times 50,000}{\text{ml sample}}$

Module No:**Module Title:**

Alkalinity Analysis

Approx. Time:

0.5 hours

Submodule Title:**Topic:**

Calculations

Instructional Objective:

Upon completion of this module the participant should be able to:

1. From titration data, determine the phenolphthalein and total alkalinity of a water sample.
2. From the data in 1., determine the hydroxide, carbonate, and bicarbonate alkalinity of a water sample.
3. Properly report the alkalinity of the water sample.

Instructional Aids:

Transparency Ak6 - Alkalinity Relationships

Instructional Approach:

Lecture/discussion

References:

Standard Methods, pp. 280-282

Class Assignments:

Module No: Ak	Topic: Calculations
Instructor Notes:	Instructor Outline:
Transparency Ak-6 Alkalinity relationships	<ol style="list-style-type: none"><li data-bbox="727 352 1584 682">1. Alkalinity: Alkalinity(CaCO_3) mg/l= $\frac{(\text{ml acid}) \times N(\text{HCl}) \times 50,000}{\text{ml sample}}$<p>ml acid is that needed to reach phenolphthalein alkalinity; and is ml acid needed to reach the total end point.</p><li data-bbox="727 741 1584 840">2. Alkalinity types: From data and table determine alkalinity of three types.<li data-bbox="727 898 1584 966">3. In reporting alkalinity, type and method should both be reported.

Exam Questions

Alkalinity Analysis

Definition of alkalinity

1. Write a chemical equation which describes how the hydroxide ion reacts with the hydrogen ion to form water.
2. Alkalinity can be used in water softening along with hardness to determine:
 - a. the rustiness of the water
 - b. the bacteria level
 - c. the relative amounts of lime and soda ash to be used
 - d. the potability of the water
3. The pH end point for phenolphthalein alkalinity is:
 - a. 8.3
 - b. 5.1
 - c. 4.5
 - d. depends on conditions
4. Which of the following ions is titrated first in an alkalinity determination?
 - a. bicarbonate
 - b. hydroxide
 - c. carbonate
5. The larger the pH value, the _____ the hydrogen ion concentration.

Use of the pH meter

6. The pH meter measures:
 - a. color
 - b. calcium ion concentration
 - c. temperature
 - d. hydrogen ion concentration
7. What solution is used to standardize the pH meter?
8. What type of electrode is placed in solution to measure pH?

Safety

9. What device in the lab can be used to wash off acid spilled all over your body?
10. _____ can be used to protect hands from acid burns.
11. Why should the cord of the pH meter be grounded?

Preparation of Reagents

12. What is the Normality of a Na_2CO_3 solution in which 2.5 g have been dissolved in water for a total volume of 1000 ml?

13. To prepare a solution which is approximately 0.1 N in HCl how many milliliters of concentrated hydrochloric acid (12 N) should be diluted to 1 liter?
- 20 ml
 - 1 ml
 - 8.3 ml
 - 100 ml

Standardization of Acid

14. When a Na_2CO_3 solution is titrated with 0.1 N HCl, what instrument is used to determine the points on the titration curve?
15. How many inflection points will be observed in the titration curve of Na_2CO_3 ?
16. If 20 ml of HCl solution are required to titrate 40 ml of a Na_2CO_3 solution which contains 2.5 g Na_2CO_3 per liter, what is the normality of the HCl solution?
17. What is the CaCO_3 equivalence of a 0.1 N HCl solution?

Titration of water samples

18. In the potentiometric titration, the first end point is due to phenolphthalein alkalinity, the second end point is due to _____ alkalinity.
19. When a water sample is titrated to a predetermined pH, what two pieces of data should be recorded?
20. What concentration acid should be used for titration of low alkalinity?

Calculations

21. If it required 10 ml of acid which is equivalent to 5 mg/ml CaCO_3 to titrate to pH 8.3, calculate the phenolphthalein alkalinity of the 50 ml sample in mg/l CaCO_3 .
22. A sample has a phenolphthalein alkalinity of zero(0) and a total alkalinity of 100 mg/l as CaCO_3 , calculate the bicarbonate alkalinity.
23. Name 3 pieces of data which should be recorded when reporting alkalinity.

ALKALINITY ANALYSIS

EQUIPMENT AND SUPPLIES LIST

1. sodium carbonate Na_2CO_3
2. drying oven
3. calcium chloride
4. dessicator
5. analytical balance
6. weighing bottle
7. 1 liter volumetric flask
8. distilled water
9. 10 ml graduated pipet
10. concentrated hydrochloric acid
11. 100 ml pipet
12. 2 - 50 ml burets
13. pH meter and electrodes
14. pH 7 buffer
15. wash bottle
16. 250 ml beaker
17. 100 ml graduated cylinder
18. bunsen burner, ring stand
19. graph paper
20. 20 ml pipet
21. magnetic stirrer and

Alkalinity Analysis

Laboratory Procedure

I. Preparation of Reagents and Standards

- A. Obtain the equipment, supplies, and chemicals listed in the "equipment" handout.
- B. Prepare the following solutions:
 1. 0.05N Na_2CO_3 . In a weighing bottle dry 5g primary standard sodium carbonate for 4 hours in an oven at 250°C . Cool in a desiccator containing CaCl_2 . Weigh the bottle on an analytical balance. Transfer 2.5g to a 1L volumetric flask and fill to the mark with distilled water. Reweigh the bottle. The mass of the Na_2CO_3 is equal to the difference of the two weighings. Call this A.
 2. 0.1N HCl. Pipet 8.3 ml concentrated hydrochloric acid (HCl) into a 1L volumetric flask. Dilute to the mark with distilled water. Mix.
 3. 0.02N HCl. Transfer 200 ml of the 0.1N HCl to a 1 liter volumetric flask. Dilute to the mark with distilled water.

II. Standardizations

A. 0.1N HCl

Fill a 50 ml buret with Na_2CO_3 solution. Fill another buret with 0.1N HCl. Standardize the pH meter with 7.0 pH buffer. Rinse the electrodes. Add 40.00 ml Na_2CO_3 solution to a 250 ml beaker. Add 60 ml distilled water. Insert pH electrodes. With constant stirring add HCl to a pH of 5.0. Remove electrodes, rinse into beaker. Boil contents of beaker for 5 minutes. Allow beaker to cool to room temperature. Titrate further, 0.2 ml at a time. Plot ml vs. pH on graph paper. Determine the inflection point (point of greatest slope). Report total ml required to reach inflection from initial minus final buret readings. Repeat procedure twice. Calculate the normality from the formula.

$$N = \frac{A \times B}{53.0 \times C}$$

Where: A is the g Na_2CO_3 used.
B is the ml Na_2CO_3 used.
C is the ml acid used.

Calculate the average N.

B. 0.02N HCl

Fill a 50 ml buret with Na_2CO_3 solution. Fill another buret with 0.02N HCl solution. Add 15 ml Na_2CO_3 solution to a 250 ml beaker. Add 65 ml distilled water. Titrate as in A above potentiometrically in triplicate. Record similar data and calculate normality (N) according to the same formula.

III. Potentiometric titration curve

- A. Preparation. Fill a 50 ml buret with standardized 0.1N HCl. Rinse the electrodes of a standardized pH meter. Pipet 20 ml well-mixed sample containing 50-200 mg total alkalinity into a 250 ml beaker. Insert electrodes.
- B. With constant stirring, add 0.2 ml increments of acid. Record the stabilized pH for each increment as well as the buret reading. Continue to pH 3.7.
- C. Prepare graph paper. Label the x-axis "ml acid added" and mark increments from 0 to the final buret reading. Label the y axis "pH" and mark in increments from 0 to 12. Title the graph "Titration curve for Alkalinity Analysis". Plot points recorded in B and connect with a smooth curve. Identify phenolphthalein inflection point (about pH 8.3) and total inflection point (about pH 5.0). Record ml HCl required to reach each end point.

IV. Titration of low alkalinity sample.

- A. Transfer 100 ml low alkalinity sample to a 250 ml beaker. Insert pH electrodes. Fill a 10 ml microburet with 0.02N standard HCl. Add acid dropwise with stirring until the pH reads 4.5.
- B. Record the exact pH and the exact number of milliliters acid required to reach this pH. Add acid dropwise with stirring to reach a pH value exactly 0.3 units less than the first pH. Record the new buret reading.

V. Calculations

- A. Potentiometric titration curve (high alkalinity)
 1. To obtain phenolphthalein alkalinity, multiply the ml acid required to reach end point by 50,000 and by the exact acid normality and divide by 20.0, the ml sample used.
 2. To obtain total alkalinity, multiply the ml acid required to reach end point by 50,000 and by the exact acid normality and divide by 20.0, the volume sample used.


3. Choose which of the following may be the case:
P=0, P less than half T. P equals $\frac{1}{2}T$, P greater than half T, P equals T. Choose the appropriate horizontal row on the data sheet and calculate hydroxide, carbonate, and bicarbonate alkalinity where appropriate. (P equals phenolphthalein alkalinity. T equals total alkalinity).
4. Calculation of low alkalinity: Double the ml titrant required to reach initial pH subtract ml required to reach final pH. Multiply this result by the acid normality and 50,000. Divide by 100, the sample volume.
5. Comment on sample source, possible errors, and other suspected ions percent.

Attach titration curves circle inflection points.

ml to inflection	I	II	III
final HCl buret reading	_____ ml	_____ ml	_____ ml
initial HCl buret reading	_____ ml	_____ ml	_____ ml
ml Acid used	_____ ml	_____ ml	_____ ml
$N = (A \times B) / (53 \times C)$	_____ N	_____ N	_____ N
Average N(HCl)	_____ N		

Titration for curve: (attach plotted graph)

ml acid	pH	ml acid	pH	ml acid	pH
---------	----	---------	----	---------	----



ml to reach end point (high alkalinity sample)

phenolphthalein _____ ml A_p
 total _____ ml A_t

Titration of low alkalinity sample.

first pH _____
 ml required _____ ml B
 second pH _____
 new buret reading _____ ml C

Calculations

High Alkalinity

phenolphthalein alkalinity
 $(A_p \times N \times 50,000) / 20.0 =$ _____ mg/l as $CaCO_3$

total alkalinity
 $(A_t \times N \times 50,000) / 20.0 =$ _____ mg/l as $CaCO_3$

Case	OH alkalinity	CO_3^{2-} alkalinity	HCO_3^- alkalinity
$P=0$	0	0	T
$P < \frac{1}{2}T$	0	2P	T-2P
$P = \frac{1}{2}T$	0	2P	0
$P > \frac{1}{2}T$	2P-T	2(T-P)	0
$P=T$	T	0	0

Bicarbonate alkalinity = _____ = _____ mg/l
 (formula) $CaCO_3$
 Carbonate alkalinity = _____ = _____ mg/l
 (formula) $CaCO_3$
 Hydroxide alkalinity = _____ = _____ mg/l
 (formula) $CaCO_3$

Low alkalinity:
 total alkalinity = $\frac{(2B-C) \times N \times 50,000}{100.0} =$ _____ mg/l
 as $CaCO_3$

Comments:

Analyst _____

Date _____

TRANSPARENCY AKI

Titration end points


total alkalinity	pH of total alkalinity end point	pH of phenolphthalein end point
30 mg/l	5.1	8.3
150 mg/l	4.8	8.3
500 mg/l	4.5	8.3
silicates or phosphates	4.5	8.3
complex systems	3.7	8.3

TRANSPARENCY AK2


Relation of pH to Hydrogen and Hydroxide Ion Concentration (Molarity)

pH	OH^-	H^+	
1	10^{-13}	10^{-1}	A C I D I C
2	10^{-12}	10^{-2}	
3	10^{-11}	10^{-3}	
4	10^{-10}	10^{-4}	
5	10^{-9}	10^{-5}	
6	10^{-8}	10^{-6}	
7	10^{-7}	10^{-7}	
8	10^{-6}	10^{-8}	B A S I C
9	10^{-5}	10^{-9}	
10	10^{-4}	10^{-10}	
11	10^{-3}	10^{-11}	
12	10^{-2}	10^{-12}	
13	10^{-1}	10^{-13}	
14	1	10^{-14}	

TRANSPARENCY AK3
Diagram of pH Meter




TRANSPARENCY AK4
Sample Na_2CO_3 Titration Curve



TRANSPARENCY AK5

Water Sample Titration Curve



TRANSPARENCY AK6

Alkalinity Relationships

P = phenolphthalein alkalinity (first inflection)

T = total alkalinity (second inflection)

	OH ⁻ alkalinity	HCO ₃ ⁻ alkalinity	CO ₃ ²⁻ alkalinity
<u>P=0</u>	<u>0</u>	<u>0</u>	<u>T</u>
<u>P 1/2T</u>	<u>0</u>	<u>2P</u>	<u>T-2P</u>
<u>P=1/2T</u>	<u>0</u>	<u>2P</u>	<u>0</u>
<u>P 1/2T</u>	<u>2P-T</u>	<u>2(T-P)</u>	<u>0</u>
<u>P=T</u>	<u>T</u>	<u>0</u>	<u>0</u>