

**Project: IEEE P802.15 Working Group for Wireless Personal Area Networks
(WPANs)**

Submission Title: [Zarlink response to 802.15 TG6 Call for Applications]

Date Submitted: [18 March, 2008]

Source: [Didier Sagan] Company [Zarlink]

Address [15822 Bernardo Center Drive, San Diego, CA 92127]

Voice:[+1 858 675-3400], FAX: :[+1 858 675-3450], E-Mail:[didier.sagan@zarlink.com]

Re: [n/a]

Abstract: [This document is Zarlink response to 802.15 TG6 Call for Applications]

Purpose: [This document is a response to P802.15 TG6 Call for Application on 18 Jan, 2008]

Notice: This document has been prepared to assist the IEEE P802.15. It is offered as a basis for discussion and is not binding on the contributing individual(s) or organization(s). The material in this document is subject to change in form and content after further study. The contributor(s) reserve(s) the right to add, amend or withdraw material contained herein.

Release: The contributor acknowledges and accepts that this contribution becomes the property of IEEE and may be made publicly available by P802.15.

Zarlink response to 802.15 TG6 Call for Applications


**Didier Sagan
Zarlink Semiconductor**

Outline

- Implanted medical devices
 - The MICS Band
 - Applications for Medical Devices
 - Design Challenges
 - *Implantable Transceiver*
- In-vitro medical devices
 - *Swallowable Camera Pill Transmitter*
- Body worn medical devices
 - *ULP Transceiver for Hearing Aids*
- Conclusion

History – Implanted Medical Telemetry

- 1980s – Inductive Telemetry
 - Near field (sub 1 MHz) at data rates <50 kHz
 - Low power (<1 mA)
 - Pick up in implant using small coil
 - Very short range (10 cm max) requiring close skin contact
- 1999 – RF Telemetry
 - Medical Implant Communication Service (MICS) Band
 - 2003 Biotronik release MICS device (non-compliant)
 - 2004 Medtronic release MICS device
 - 2005 Guidant release ISM band (915 MHz) device
 - 2007 St-Judes release MICS device
 - ISM bands (13.56, 433, 868, 915 MHz) are sometimes used
- 2002 - Ultrasonic Telemetry




What is MICS ?

- Medical Implant Communication Service (MICS)
 - 402–405 MHz frequency allocation
 - FCC was petitioned in 1988, allocated in 1999
 - Short-range, wireless link to connect low-power implanted medical devices with monitoring and control equipment
 - Implanted Medical Devices (IMD) such as cardiac pacemakers, implantable cardioverter/defibrillator (ICD), neurostimulators, etc.
 - Why 402-405 MHz?
 - Reasonable signal propagation characteristics in the human body
 - Compatibility with incumbent users of the band (e.g. weather balloons)
 - General world-wide acceptance
 - Approved in United States, Europe, Canada, Australia and Japan

Why was MICS Introduced?



- Traditional telemetry not user friendly
 - Use inductive links
 - In contact with patient
 - Very limited range
 - Low frequency
 - Data rates similar to a dial-up computer modem
 - Not user friendly for home monitoring
 - Requires a wand to be positioned above the IMD by the patient
- Need for higher data rates
 - To upload patient events captured in the IMD's memory to the base station for analysis
 - Reprogram the IMD
 - Shorten doctor/patient consultancy times



Why was MICS Introduced?

- Need for longer range
 - Locate the base station (programmer) outside of the sterile field during surgery
 - Simplify home-monitoring for elderly
 - Broaden possible applications
 - Bedside monitor for emergency


- Competitive pressures of the medical device industry
 - Higher data rates enable new, value-added services



MICS - Applications

- **Stimulatory Devices**
 - Pacemaker
 - Implantable Cardioverter/Defibrillator (ICD)
 - Neurostimulators and pain suppression devices
 - FES (Functional Electrical Stimulation)

- **Measurement/Control/Other Devices**
 - Drug infusion and dispensing
 - Artificial heart and heart assist devices
 - Implanted sensors
 - Control of other artificial organs and implanted devices




Challenges

- Ultra Low Power Consumption
 - Facts about implanted pacemakers
 - Lifetime > 7 years; up to 10 years
 - Maximum current drain of the order of 10 – 20 μA
 - Telemetry budgeted as no more than 15%, i.e. 2 – 3 μA
 - Telemetry is off most of the time but still need to sniff every 1 – 10 s
 - Consumption during Sleep/Sniff modes is therefore the most critical
 - Requirements
 - Very low TX/RX current, <15nA.s/bit
 - Ultra low sleep/listen current, ideally <100s of nA
 - Implies leakage <100nA (at body temperature)



Challenges

- **Minimum External Components**
 - RF module <3x5x10 mm
 - Fewer components => higher reliability, lower cost, smaller size
- **High data rates**
 - Need to go beyond current payload requirements
 - Sending data in short bursts conserves power
 - Reduces time window for interference and easier supply decoupling
- **Operating range**
 - Require ~2 m outside the body to improve on existing links (short range inductive)
 - Antenna matching, fading and body loss typically 40-50 dB
- **Reliability**
 - Data and link integrity, selectivity and interference rejection
 - Need to work 24/7 for 10 years; no reset button...



Module size 3 x 5 x 10 mm


ULP Implantable Transceiver

MICS and ISM Band Transceiver:

- Negligible standby current
- high data and low error rates in a small footprint



ZL70101:

Supply Voltage:	2.1 - 3.5 V Battery
Radio Frequency:	402-405 MHz / 433 MHz
Type of RF link	Bi-directional, half duplex
Modulation Scheme:	FSK
Raw Bit Rate:	800 / 400 / 200 kbits/s
Operating Current:	5mA TX/RX down to <1mA
Sleep (sniff) Current:	< 250 nA
Ext. comps:	3 (2 caps, Xtal) + antenna matching)
BER:	< 1.5×10^{-10}
Range:	~2 m
Interface:	SPI



Wakeup

- Problem: MICS band limited to 25 μW (-16 dBm)
- Consequence: Very small received signal
 - Receiver power too large to meet both latency and power consumption requirement in sniff mode.
- One solution: Use band with more power 2.45 GHz (up to 20 dBm) and design a synthesizer-less receiver
 - 250 nA average current for 1.15 second latency



The Camera Pill

(Company: Given Imaging)

New digestive track diagnostic tool


Replace traditional endoscopy

Better diagnostic

3D mapping capability



Healthy Small Bowel



The Camera Pill - Facts



Size: 11 x 26 mm

Weight: < 4 gram


View: 140 deg

Approximately 57,000 pictures

During 8 hours

The Camera Pill – The Inside

- World's First Swallowable Camera Capsule, from Given Imaging, including Zarlink's ULP RF Transmitter




A real "Fantastic Voyager," Given's tiny camera capsule is swallowed.

ULP Medical Transmitter

Very high data rate transmitter

- ➔ very low power
- ➔ small footprint
- ➔ designed for imaging applications

Supply Voltage	2.6 - 3.2 V Battery
Radio Frequency:	400 - 440 MHz
Type of RF link:	Transmit only
Bit Rate:	2700 kbits/s
Operating Power:	5.2 mW
Ext. comps:	10



Hearing Aids are becoming Communication Devices

- Programming would be easier, more reliable and cheaper without cable.
- Ear-to-ear communication can improve hearing by coming closer to providing a real stereo image.
- Streaming audio allow connection to cell phone and MP3 player.

- But power consumption is limited
 - HA need to last 1 – 2 weeks, 16h/day
 - Battery capacity is limited (typ ~250 mA.h, 1.25V)
 - Battery chemistry limits peak current (few mA)
 - HA functionality is ~1mA
 - Wireless link budget should be < 25%



Physiology Monitoring

- Many different use cases
 - Clinical environment
 - Home monitoring (chronical disease, health status)
 - Fitness
 - Animals (Lab test, pets & live stock)
- Parameters:
 - HR, Temp, BP, ECG, EEG, SpO2, Plethysmography, etc.
- Usage model:
 - Spot measurements vs Continuous monitoring
 - From few days to few months to few years
- Number of sensors: 1 to 6
- Payload:
 - From few bits/s up to 5 kbps per sensor
 - Aggregate: up to 10-15 kbps
- In some cases a wireless node would need ~10kbps, operating continuously 24/7 for months without changing the battery.
 - Wireless telemetry needs an efficiency better than 15nJ/bit

ULP ISM Transceiver

Requirements:

- Smallest possible peak power
- Data rate that supports audio streaming
- Very efficient (11 nJ/bit in both Tx and Rx)
- Very small



ZL70250:

Radio Frequency: 915 MHz (Americas) / 863-865 MHz (Europe)

Type of RF link: Bi-directional, half duplex

Bit Rate: 186 kbits/s

Current Consumption: 2 mA from 1.05 - 1.5 V Battery

Range: 3 meters (limited by antenna size)

Externals: 2 (Xtal, Res)

Interface: SPI + 2-wire

Comments about antennas

- Assumption that antenna size shrinks with higher frequency is not always true
 - For given range and path loss, antenna size remains about the same
 - In other materials than air, the practical antenna type vs freq combinations can compensate or even reverse the relation.
- The dipole is not the best antenna near the body
 - A dipole, like any electro-magnetic antenna, is detuned by the body
 - An electrically small loop is behave as a magnetic antenna and is very little affected by the body
 - But electrically small loops are not practical for GHz frequencies
- Practical example:
 - 25mm x 10mm loop at 900MHz
 - Built on FR4 PCB: very cheap and repeatable
 - -5dBi antenna gain
 - Not affected by pinching it while same (crude) test with a dipole resulted in significant radiation loss.

Conclusion

- RF telemetry in and around the body already exists. There are commercial parts available to support such applications
- For implanted telemetry, most of the industry has chosen to use the MICS band.
- Most 'around the body' applications use a proprietary solutions, often by choice (seen as a competitive differentiator)
- Enabling parameter: ULP power consumption

